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Development and Validation of Simultaneous Determination Method of Benfuresate and Cyhalofop-butyl in Whole-crop Rice Silage by LC-MS/MS

Kozue SATOU^{*1}, Eiko ARAI^{*1} and Reiko KAZAMA^{*2}

(^{*1} Sendai Regional Center, Food and Agricultural Materials Inspection Center (FAMIC)

^{*2} Fertilizer and Feed Inspection Department, FAMIC)

An analytical method was developed to simultaneously determine the levels of benfuresate and cyhalofop-butyl in whole-crop rice silage (WCRS) using a liquid chromatograph-atmospheric pressure chemical ionization-tandem mass spectrometer (LC-APCI-MS/MS).

After adding water to the sample, benfuresate and cyhalofop-butyl were extracted with acetone and the extracted solutions were filtered. The filtrate was purified with two types of SPE columns (InertSep K-solute, GL Sciences Inc.; Tokyo, Japan and ENVI-Carb/LC-NH₂, Sigma-Aldrich Co. LLC.; St. Louis, MO, USA), and injected into the LC-APCI-MS/MS to determine the levels of benfuresate and cyhalofop-butyl. LC separation was carried out on an ODS column (Inertsil ODS-3, 2.1 mm i.d. × 150 mm, 3 μm from GL Sciences Inc.) using a gradient with 0.1 v/v% formic acid aqueous solution and 0.1 v/v% formic acid methanol solution as a mobile phase. In the MS/MS analysis, positive mode atmospheric pressure chemical ionization (APCI+) was used.

Recovery tests were conducted on WCRS spiked with benfuresate at the levels of 0.009 or 0.2 mg/kg original matter and cyhalofop-butyl at the levels of 0.009 or 0.1 mg/kg original matter. The resulting mean recoveries ranged from 86.0 % to 101 % for benfuresate and 85.3 % to 93.2 % for cyhalofop-butyl. The repeatability in terms of the relative standard deviations (RSD_r) were not more than 12 % for benfuresate and 14 % for cyhalofop-butyl.

A collaborative study was conducted in eleven laboratories using two kinds of WCRS. One WCRS was spiked with 0.2 mg/kg air-dry matter of benfuresate and 0.1 mg/kg air-dry matter of cyhalofop-butyl. The other of WCRS was spiked with 0.02 mg/kg air-dry matter of benfuresate and cyhalofop-butyl. The resulting mean recoveries, RSD_r, reproducibility in terms of the relative standard deviations (RSD_R), and HorRat, respectively for two kinds of spiked WCRS, were 96.1 % and 94.7 %, 2.6 % and 3.7 %, 8.1 % and 10 %, and not more than 0.45 for benfuresate, and 92.7 % and 97.0 %, 2.7 % and 4.8 %, 22 % and 32 %, and not more than 1.5 for cyhalofop-butyl.

This method was validated and established for use in the inspection of benfuresate and cyhalofop-butyl in whole-crop rice silage.

Key words: benfuresate; cyhalofop-butyl; liquid chromatograph-tandem mass spectrometer (LC-MS/MS); atmospheric pressure chemical ionization (APCI); whole-crop rice silage; collaborative study

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

Benfuresate

2,3-dihydro-3,3-dimethylbenzofuran-5-yl ethanesulfonate

$C_{12}H_{16}O_4S$ MW: 256.3 CAS No.: 68505-69-1

Cyhalofop-butyl

butyl (*R*)-2-[4-(4-cyano-2-fluorophenoxy)phenoxy]propionate

$C_{20}H_{20}FNO_4$ MW: 357.1 CAS No.: 122008-85-9

Fig. 1 Chemical structures of benfuresate and cyhalofop-butyl

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
Fig. 4-1 Calibration curves of benfuresate by peak area (left) and peak height (right)


Fig. 4-2 Calibration curves of cyhalofop-butyl by peak area (left) and peak height (right)

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Table 3 Effects of volume of sample solution for purification on recoveries

Pesticides	Spiked level (mg/kg air-dry matter)	Recovery (%)			
		Volume of sample solution for purification			
		5 mL	10 mL	20 mL	40 mL
Benfuresate	0.4	95.8	116	120	116
Cyhalofop-butyl	0.2	98.3	106	86.3	44.6

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Table 4 Elution pattern of benfuresate and cyhalofop-butyl from InertSep K-solute

Pesticides	Recovery (%)				Total
	Hexane				
	0~60 mL	60~80 mL	80~100 mL	100~120 mL	
Benfuresate	90	0	0	0	90
Cyhalofop-butyl	90	2	1	0	93

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Pesticides	Recovery (%)			Total
	Ethyl acetate			
	0~10 mL	10~12 mL	12~14 mL	
Benfuresate	99	0	0	99
Cyhalofop-butyl	100	0	0	100

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
Fig. 5 Selected reaction monitoring (SRM) chromatograms of standard and blank sample solution

(LC-MS/MS conditions are shown in Tables 1 and 2. Arrows indicate the retention times of 1: benfuresate and 2: cyhalofop-butyl. The baselines are shifted for display.)

A: Standard solution (5 ng/mL: 50 pg as each pesticide)

B: Sample solution of WCRS (blank)

Table 6 Results of matrix effect study

Pesticides	Concentration of pesticides		Matrix effect ^{b)} (%)
	in matrix standard solution (ng/mL)	in sample ^{a)} (mg/kg air-dry matter)	
Benfuresate	5	0.02	100
	100	0.4	100
Cyhalofop-butyl	5	0.02	112
	50	0.2	107

$n = 1$

a) Converted from the concentration in matrix standard solution

b) Ratio of peak area of pesticides in the presence of matrix to that in the absence of matrix

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a) Mean ($n = 5$)

b) Relative standard deviation of repeatability

c) The pesticides were spiked to air-dried WCRS samples one night prior to extraction. The spiked levels were 0.02 and 0.4 mg/kg air-dry matter for benfuresate and 0.02 and 0.2 mg/kg air-dry matter for cyhalofop-butyl, respectively. The levels of pesticides in original matter were calculated with following equation on the assumption that the moisture content of WCRS samples was 60 % for original matter and 10 % for air-dry matter.

The levels of pesticides in original matter (moisture 60 %)

= the levels of pesticides in air-dry matter (moisture 10 %) / 2.25


Fig. 6 SRM chromatograms on recovery test

(LC-MS/MS conditions are shown in Tables 1 and 2. Arrows indicate the peaks of 1: benfuresate and 2: cyhalofop-butyl. The baselines are shifted for display.)

A: Standard solution (100 ng/mL: 1 ng as benfuresate and 50 ng/mL: 0.5 ng as cyhalofop-butyl)

B: Sample solution of WCRS (spiked at 0.2 mg/kg original matter of benfuresate (as 100 ng/mL in sample solution), and 0.1 mg/kg original matter of cyhalofop-butyl (as 50 ng/mL in sample solution))

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Table 8 Collaborative study for benfuresate

Lab. No.	Sample			
	WCRS 1		WCRS 2	
	(mg/kg air-dry matter)		(mg/kg air-dry matter)	
1	0.426	0.420	0.0435	0.0407
2	0.398	0.397	0.0386	0.0394
3	0.459	0.437	0.0466	0.0512
4	0.444	0.408	0.0470	0.0456
5	0.403	0.390	0.0394	0.0375
6	0.480	0.483	0.0370	0.0379
7	0.442	0.445	0.0413	0.0409
8	0.485	0.509	0.0475	0.0482
9	0.383	0.394	0.0402	0.0375
10	0.444	0.435	0.0468	0.0458
11	0.420	0.415	0.0440	0.0409
Spiked level (mg/kg air-dry matter)	0.45		0.045	
Mean value ^{a)} (mg/kg air-dry matter)	0.433		0.0426	
Mean recovery ^{a)} (%)	96.1		94.7	
RSD _r ^{b)} (%)	2.6		3.7	
RSD _R ^{c)} (%)	8.1		10	
PRSD _R ^{d)} (%)	18		22	
HorRat	0.45		0.45	

a) $n = 22$

b) Relative standard deviation of repeatability within laboratory

c) Relative standard deviation of reproducibility between laboratories

d) Predicted relative standard deviation of reproducibility between laboratories calculated from the modified Horwitz equation

Table 9 Collaborative study for cyhalofop-butyl

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a) $n = 22$

b) Relative standard deviation of repeatability within laboratory

c) Relative standard deviation of reproducibility between laboratories

d) Predicted relative standard deviation of reproducibility between laboratories calculated from the modified Horwitz equation

Table 10 Instruments used in the collaborative study

Lab.No	LC-MS/MS	LC column (i.d.×length, particle size)
1	LC: ACQUITY UPLC, Waters	Inertsil ODS-3, GL Sciences
	MS/MS: Quattro premier XE, Waters	(2.1 mm×150 mm, 3 μm)
2	LC: Nexera X2, Shimadzu	Inertsil ODS-3, GL Sciences
	MS/MS: LCMS-8040, Shimadzu	(2.1 mm×150 mm, 3 μm)
3	LC: Nexera X2, Shimadzu	Inertsil ODS-3, GL Sciences
	MS/MS: LCMS-8040, Shimadzu	(2.1 mm×150 mm, 3 μm)
4	LC: ACQUITY UPLC, Waters	Inertsil ODS-3, GL Sciences
	MS/MS: ACQUITY TQD, Waters	(2.1 mm×150 mm, 3 μm)
5	LC: ACQUITY UPLC, Waters	ZORBAX Eclipse XDB-C18
	MS/MS: ACQUITY TQD, Waters	Agilent Technologies (4.6 mm×150 mm, 3.5 μm)
6	LC: ACQUITY UPLC, Waters	Mightysil RP-18GP
	MS/MS: ACQUITY TQD, Waters	Kanto Chemical (2.0 mm×150 mm, 3 μm)
7	LC: Nexera XR, Shimadzu	Inertsil ODS-3, GL Sciences
	MS/MS: QTRAP 5500 AB SCIEX	(2.1 mm×150 mm, 3 μm)
8	LC: Alliance 2695, Waters	Inertsil ODS-3, GL Sciences
	MS/MS: Quattro premier XE, Waters	(2.1 mm×150 mm, 3 μm)
9	LC: Prominence Shimadzu	Inertsil ODS-3, GL Sciences
	MS/MS: QTRAP 5500 AB SCIEX	(2.1 mm×150 mm, 3 μm)
10	LC: Nexera X2, Shimadzu	Inertsil ODS-3 HP, GL Sciences
	MS/MS: LCMS-8040, Shimadzu	(2.1 mm×150 mm, 3 μm)
11	LC: 1200 Series Agilent Technologies	Inertsil ODS-3, GL Sciences
	MS/MS: 6410 Triple Quadrupole LC/MS Agilent Technologies	(2.1 mm×150 mm, 3 μm)

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- 8) William Horwitz: Protocol for design, conduct and interpretation of method-performance studies, *Pure & Appl. Chem.*, **67**(2), 331-343 (1995).
- 9) George W. Latimer, Jr.: Official methods of analysis of AOAC INTERNATIONAL 20th edition, Appendix D, Guidelines for collaborative study procedures to validate characteristics of a method of analysis. Gaithersburg, MD, USA (2016) (ISBN: 978-0-935584-87-5).
- 10) Michael Thompson: Recent trends in inter-laboratory precision at ppb and sub-ppb concentrations in relation to fitness for purpose criteria proficiency testing, *Analyst*, **125**, 385-386 (2000).